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which was used in laboratories. The results below were obtained by 159 subjects with diabetes at three independent clinical sites. The regression statistics are derived from a plot of the smartLAB® nG capillary data versus YSI plasma data

Slope	0.9803
Intercept	+ 2.9625 mg/dL (+ 0.1646 mmol/L)
Correlation coefficient	13.6852
Number of subject	159
Range tested	47.5–566.5 mg/dL or 2.64–31.74 mmol/L
R-square	0.9782

REFERENCES  
1. *Cass, A.E.G. et al., Anal. Chem., 56 (1984) p.667*  
2. *Tietz N.: Fundamentals of clinical chemistry, 3rd Ed., W.B. Saunders Co., Philadelphia, PA, 1987, p.427*  
3. *NCCLS document M29-A, Protection of the laboratory worker from instrument biohazards and infectious disease transmitted by blood, body fluids, and tissue, 1997.*  
4. *American Diabetes Association. Diagnosis and Classification of Diabetes Mellitus (Position Statement). Diabetes Care 34 (Suppl. 1) S66, 2011.*  
5. *American Diabetes Association. Standards of medical care in diabetes. Diabetes care. 2013; Vol. 36, Suppl 1, S21*

- LIST OF SYMBOLS
-  Caution. Please read manual thoroughly.  In-vitro diagnostics only. Not to ingest.
  -  Expiration date  Storage temperature range
  -  This product meets the requirements of directive 98/79 CE for in-vitro diagnostic medical devices.  Control solution
  -  Please read instructions thoroughly.  Opening date of test strip via.
  -  Sufficient for  Single use only.
  -  Manufactured by  Lot number
  -  Blood glucose test result in mg/dL  Catalogue number
  -  Blood glucose test result in mmol/L  6 months until the product expires upon opening

NOTE ON **smartLAB® nG** -TEST STRIPS  
The **smartLAB® nG** blood glucose test strips compatible with **smartLAB® nG** glucose meters need not to be coded, any more.

 **HMM Diagnostics GmbH** D/N D11450-08-02-20\_R2\_12/12/2020  
Schulze-Delitzsch-Str. 16   
D-68542 Heddeshheim, Germany  
www.hmm.info 

Blood Glucose System.  
5. Interfering Substances depend on the concentration. The below substances up to the test concentration will not affect the readings.

Concentrations of the interference tested	Bias		Glucose Level	80 mg/dL (4.4 mmol/L)	250 mg/dL (13.9 mmol/L)	500 mg/dL (27.8 mmol/L)
Ascorbic Acid	4 mg/dL	(0.26 mmol/L)		10.89%	-1.76%	4.55%
Ibuprofen	50 mg/dL	(2.43 mmol/L)		3.10%	2.88%	4.62%
L-Dopa	1.8 mg/dL	(0.09 mmol/L)		10.59%	7.91%	4.90%
Sodium Salicylate	50 mg/dL	(3.12 mmol/L)		-2.59%	9.42%	-0.84%
Tetracycline	1.5 mg/dL	(0.03 mmol/L)		-5.32%	3.81%	3.20%
Tolbutamide	100 mg/dL	(3.70 mmol/L)		-2.60%	12.30%	0.89%
Bilirubin-unconjugated	2.4 mg/dL	(0.04 mmol/L)		-2.52%	4.05%	-0.23%
Uric acid	8 mg/dL	(0.48 mmol/L)		2.71%	9.55%	-1.75%
Xylose	4 mg/dL	(0.27 mmol/L)		-5.12%	-1.64%	-4.44%

Use the anticoagulant EDTA or sodium Heparin or Lithium Heparin to preserve the venous blood and capillary blood specimens would not affect the readings of **smartLAB® nG** System. If the specimen was preserved by the anticoagulant NaF/potassium Oxalate contained sodium fluoride, the reading of **smartLAB® nG** System would be seriously interfered.

PERFORMANCE CHARACTERISTICS  
The test range is between 20 to 630 mg/dL (1.1~35.0 mmol/L). Validate the test strips performance in both laboratory and clinical tests.

PRECISION  
Precision of the **smartLAB® nG** Blood Glucose System was measured with both venous blood and control solution in the laboratory. The results are shown in the table below:

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ACCURACY  
The accuracy of the **smartLAB® nG** Self-Monitoring Blood Glucose System was assessed by comparison

questions for the readings, check the following items first and then repeat the test. If the results are still questionable, consult your doctor:  
1. If the strips are within the expiration date.  
2. Make sure the drop of blood in the whole reaction zone.  
3. Check meter and test strip performance with **smartLAB® nG** Glucose Control Solutions.

Expected Results:  
The **smartLAB® nG** Blood Glucose test strips are plasma referenced and calibrated for easier comparison for lab results.  
• For Non-Diabetic\*  
The American Diabetes Association recommends a post-meal glucose level of less than 140 mg/dL (7.8 mmol/L) and a pre-meal glucose of less than 100 mg/dL (5.6 mmol/L).  
• For Non-Pregnant Adults with diabetes<sup>5</sup>  
The American Diabetes Association recommends a post-meal glucose level of less than 180 mg/dL (10.0 mmol/L) and a pre-meal glucose of 70-130 mg/dL (3.9-7.2 mmol/L).  
**CAUTION:** Any low or high blood glucose readings can indicate a potentially serious medical condition. If the readings do not reflect your symptoms, repeat the test with a new test strip. Contact your doctor when your reading is:  
A. Not consistent with your symptoms.  
B. Less than 60 mg/dL (3.3 mmol/L).  
C. Higher than 240 mg/dL (13.3 mmol/L).

LIMITATIONS  
The test strips are used for venous blood or fresh capillary whole blood samples.  
1. DO NOT use serum or plasma sample.  
2. DO NOT use anticoagulant NaF or potassium oxalate for venous sample preparation.  
3. DO NOT use neonate blood sample.  
4. Extreme humidity may affect the results. A relative humidity greater than 90% may cause incorrect results.  
5. Use this system at temperatures between 10°C and 40°C (50°F and 104°F). Outside this range, the system may get incorrect results.  
6. DO NOT reuse the test strips. The test strips are for single use only.  
7. Hematocrit (HCT): HCT between 35% - 50% will not significantly affect the results. HCT at 30% - 35% and at 50% - 55% may affect the reading with ±15%. HCT below 35% may cause higher results. HCT above 50% may cause lower results. If you do not know your HCT level, please consult with your doctor.  
8. Altitude up to 3,048 meters above sea level has no effect on readings.  
Healthcare providers – Please note these additional Limitations.

1. If the patient has the following conditions, the result may be inaccurate:  
• Severely dehydrated.  
• Severely hypotensive.(low blood pressure)  
• In shock.  
• In hypoglycemic-hyperosmolar state.(with or without ketosis)  
2. Lipemic samples: Cholesterol level up to 500 mg/dL (12.92 mmol/L) and triglycerides up to 3,000 mg/dL (33.6 mmol/L) do not affect the results. Do not test grossly lipemic patient samples with **smartLAB® nG** System.  
3. Do not use home-use blood glucose systems to test critically ill patients.  
4. DO NOT use during xylose absorption testing. Xylose in the blood will interfere Self-Monitoring

3. Consult with your healthcare professional before you test palm and forearm.  
TESTING YOUR BLOOD GLUCOSE LEVEL  
Wash and dry your hands thoroughly. Insert the lancet into the lancing device. Prepare the lancing device according to the user's manual. To test your blood glucose level, please follow these steps to test:

Step 1: Insert Test Strip. Take a test strip and replace the bottle cap quickly. Insert the test strip into the slot of the meter. The meter will turn on.  
Step 2: Apply Sample. Follow the sample collection procedures to get a drop of blood. While the “” appears, bring your blood to the collection area on the test strip. The blood will automatically be drawn into the reaction zone. The meter will begin to run.  
Step 3: Accurate results after 5 seconds. Your result will appear on the screen after 5 seconds.  
Step 4: Thoroughly wash hands with soap and water after performing your glucose testing. For detailed information on the test procedure, please refer to the User's Manual.

OPTIONAL: QUALITY CONTROL  
It is recommended to test the performance of the **smartLAB® nG** self-monitoring Blood Glucose System when  
1.Begin using a new bottle of test strips.  
2.Question your blood test results.  
3.Check that the meter and test strips are working properly, and that you perform the test correctly at least once.  
4.Suspect the meter or test strips are working incorrectly.  
5.Drop the meter.  
6.Expose the test strips to temperatures outside the storage conditions.(2°C to 30°C,36°F to 86°F)  
Use the control solution to check the performance of the meter and test strips. Please follow the instructions described in the section “Control Solution Test” of the User's Manual. Use the **smartLAB® nG** Glucose Control Solution to test the system. Don't use other manufacturer's control solution. Other control solution may provide incorrect results.

IMPORTANT  
If the control solution test results are outside the range, your meter and test strip may not be working as a system. Follow the correct procedures and repeat the control solution test. Make sure the control solution is clean without contamination. (e.g. Turbid.) Do not use the meter until test results fall within the appropriate range. If the problem continues, please contact your local dealer.  
When the results fall within the specified ranges printed on the strip bottle label, the system is working properly. If your control solution test results are outside the range, repeat the control solution test. Results that fall outside the range may be caused by:  
1. Expired or contaminated control solution.  
2. Error in performing the test.  
3. Meter malfunction.  
4. Test strip deterioration.  
DO NOT use the system to test your blood until you get a control solution test result within the specified range.

TEST RESULTS  
If “LO” appears on the screen, your blood glucose level is lower than 20 mg/dL (1.1 mmol/L). If “HI” appears, that your blood glucose level is higher than 630 mg/dL (35.0 mmol/L). When you get any

3. Do not re-use test strips.  
4. Do not cut, bend, scratch, or alter the test strip in any way.  
5. Warning for potential biohazard<sup>3</sup>:  
This warning is for the healthcare providers who use this system on multiple patients. They should know the product comes in contact with human blood. After cleaning, dispose of the product properly as if avoid transmitting viral disease.  
6. If the readings are not consistent to your symptoms, check the manual first. Make sure you have followed all instructions, and then contact your doctor.  
7. Never make any changes to your diabetes control program without consulting your doctor.  
8. The system should be used at a temperature between 10°C and 40°C (50°F and 104°F). Outside this range, the system may get incorrect results.

SUPPLIES FOR TESTING YOUR BLOOD GLUCOSE  
**smartLAB® nG** Blood Glucose System.  
**smartLAB® nG** Blood Glucose Test Strip.  
Lancing Device  
Lancets.  
**smartLAB® nG** User's Manual.

SAMPLE COLLECTION AND PREPARATION  
This system is for use in fresh capillary or venous whole blood. The blood sample must be used right away. The sample volume requires at least 0.6 micro-liters. To obtain a drop of blood, follow these steps:

Step 1: Wash and dry your hands thoroughly.  
Step 2: Prepare the meter according to the user's manual.  
Step 3: If you use an alcohol swab. Make sure that your finger is entirely dry before lancing.  
Step 4: Puncture and get a drop of blood. Avoid pressing too hard against the punctured site.

ALTERNATE SITE TESTING (AST)  
1.AST results may differ from fingertip results when glucose levels are changing rapidly (e.g., after a meal, after taking insulin, or during or after exercise).  
2.AST can only be used during steady-state blood glucose conditions. AST should only be performed under the following conditions:  
• Testing before a meal.  
• In a fasting state.  
• 2 hours or more after a meal.  
• 2 hours or more after insulin dosing.  
• 2 hours after physical activity.

Fingertip test only:  
• If sick.  
• If blood glucose is low.  
• After exercising.  
• When you have just taken insulin  
• Two hours or less after eating.  
• If you often do not notice when your blood glucose is low.  
• After injecting rapid-acting insulin (2 hours or less).

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3. NCCLS document M29-A, *Protection of the laboratory worker from instrument biohazards and infectious disease transmitted by blood, body fluids, and tissue*, 1997.

4. American Diabetes Association. *Diagnosis and Classification of Diabetes Mellitus (Position Statement)*. *Diabetes Care* 34 (Supp. 1) S66, 2011.

5. American Diabetes Association. *Standards of medical care in diabetes*. *Diabetes care*. 2013; Vol. 36, Suppl 1, S21

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www.hmm.info

**NO CODE**

## INSTRUCTION **smartLAB<sup>®nG</sup>** TEST STRIPS

### INTENDED USE

The **smartLAB<sup>®nG</sup>** blood glucose systems are intended for the quantitative measurement of glucose in venous blood or fresh capillary whole blood samples drawn from the fingertip, palm or forearm. Testing is done outside the body (In Vitro diagnostic use). It is indicated for self-testing by persons with diabetes, or in clinical settings by healthcare providers, as an aid to monitor the effectiveness of diabetes control. The system is not to be used on neonates, is not for the diagnosis of or screening for diabetes mellitus, and that alternate site testing can only be used during steady-state blood glucose conditions.

### IMPORTANT

Please read this insert and the **smartLAB<sup>®nG</sup>** User's Manual before using **smartLAB<sup>®nG</sup>** Blood Glucose Test Strips. If you have any questions and/or need assistance, please contact our authorized distributors in your country.

### MEASUREMENT PRINCIPLE

Apply a blood sample to the blood collection area at the tip of the strip. The test strip will draw the blood sample into the reaction zone. The FAD-binding glucose dehydrogenase catalyzes the glucose dehydrogenation. It will produce gluconolactone. During the reaction, a mediator transfers electrons to the electrode surface and generates a current. The amount of the current is proportional to the amount of glucose present in the blood sample. After 5 seconds, the **smartLAB<sup>®nG</sup>** system will show glucose concentration on the screen.

### REAGENT COMPOSITION

Each test strip contains the following reagent:

6 (w/w) %	FAD glucose dehydrogenase (Aspergillus oryzae, 2.0 IU/test strip)
56 (w/w) %	Potassium ferricyanide
38 (w/w) %	Non-reactive ingredients

### STORAGE AND HANDLING

- Store the test strip package in a cool, dry place, between 2°C~30°C (36°F - 86°F). Do not freeze.
- Keep out of direct sunlight and heat.
- Write the date on the bottle when first opened.
- Discard test strips 6 months after first opening date.
- Replace the cap quickly after taking the test strip from the bottle.
- Handle the test strip only with dry and clean hands.
- Store your test strips in their original bottle only. Do not move them to any other bottle or container

### WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only. Do not swallow.
- Do not use the test strip after the expiration date.

**smartLAB<sup>®nG</sup>**