smartLAB[®] fit

Körperanalysewaage

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Body Analyser Scale

User Manual

Please read the user manual carefully before the first use.

*Current information about this device, visit smartLAB on
www.support.hmm.info*

HMM Diagnostics GmbH
Friedrichstraße 89
D-69221 Dossenheim, Germany
www.hmm.info

Content

I. Introduction	22
II. Your smartLAB® fit	23
Display & functions	23
Specifications	24
Set Content	24
General Information	24
Security information	26
III. Setup & Operating functions	27
Changing the batteries	27
Beginning of operation / „quick weighing “	27
Setting up of the user data	27
Measuring	29
Optional: Weighing on a carpet	30
Interpreting the results	30
IV. Miscellaneous	34
Maintenance	34
Regulations and Safety Notice	34
Warranty	35

I. Introduction

Thank you for having decided in favor of our **smartLAB® fit** body analyser scale. This product enables you to measure your body weight as well as content of adipose, content of water, muscle and bone mass. Please read the instructions carefully before you use the scale for the first time and keep the instructions for a future use.

Principle of measurement


Your **smartLAB® fit** body analyser scale is based on the principle of the bioelectric impedance analysis (BIA). A weak and non-perceptible and, therefore, not dangerous electric signal will be conducted through the whole body. The scale identifies the electric resistance (impedance). Individual data such as age, body height, gender or degree of activity allow to identify the adipose content and other parameters.

Muscle tissue and water have a good electric conductivity and, therefore, a lower resistance than bone- and adipose tissue because adipose cells and bones hardly conduct current.

Please note that the values of our **smartLAB® fit** body analyser scale are only an approach to the medical, real values of the body. Only a specialist who takes into account medical methods (e. g. computer tomography) can identify the exact parameters like adipose, content of water, content of muscles and bone structure.

II. Your smartLAB® fit

Display & functions

The settings can be adjusted by pressing on the side of the scale.

Specifications

1. Device: **smartLAB[®] fit**
2. Max. capacity: 150 kg / 330 lb
3. Units: kg / lb (in 100 g / 0,2 lb steps)
4. Buttons: Set, Go up, Go down
5. Activity indices: 4
6. Adipose: in % for 0,1% steps
7. Content of water of the body: in % for 0,1 % steps
8. Muscle mass: in % for 0,1% steps
9. Basal metabolic syndrome: in kcal with 1 kcal per step
10. Active metabolic syndrome: in kcal with 1 kcal per step
11. User profiles: 8
12. Battery: 4 x AAA batteries (1,5 V)
13. Size (L x W x H): 310 x 310 x 19
14. Weight: 1950 g
15. Operating temperature: 0°C - 40°C
16. Storage temperature: -20°C - 60°C
17. Symbol when overburdened: „----“
18. Symbol for low batteries: „Lo“

Contents

- 1 **smartLAB[®] fit** diagnostic scale
- 1 User manual
- 4 1,5V „AAA“ batteries
- 4 carpet extension feet

General Information

- If possible please use the scale at the same time every day (we recommend in the morning), after having used the toilet, with empty stomach and without clothes in order to be able to compare the results.
- The adipose can only be detected bare-footed. You can moisten the soles of your feet. Completely dry feet can result in incorrect measurements because they do not have enough conductivity.
- Please do not move during the measurement.
- Please wait several hours before measuring after unaccustomed physical exercise.
- Please wait about 15 minutes after getting up before you measure your body weight so that the water of your body can disperse.
- The following persons or groups of persons will get differing or not plausible results concerning their adipose and the other values:

- children under 10 years
 - competitive athletes & bodybuilders
 - pregnant woman
 - people having fever, dialysis patients, people with edema symptoms or osteoporosis
 - persons taking cardiovascular medicine
 - persons taking vasodilatory or vasoconstrictive medicine
 - persons with considerable anatomic differences of their legs that can influence the height of the body (different lengths of the legs)
-
- The device is for your own use only and not for medical or commercial use.
 - Please note that technical related measuring tolerances are possible since the scale is not especially made for professional medical use.
 - Please place the scale on an even, hard floor. The floor is the basis for a correct measurement.
 - The **smartLAB[®] fit** body analyser scale corresponds to the EC directive 2004/108.

Please feel free to contact the manufacturer or the customer service for more information.

Security information

The **smartLAB**[®] *fit* body analyser scale must not be used by persons with medical implants (e. g. cardiac pacemaker). Their function could be affected.

Do not use during pregnancy. The amniotic liquor inside the body can influence the measurements.

Please do not use the **smartLAB**[®] *fit* body analyser scale with wet feet – slip hazard!

Do not swallow the batteries – danger of life. Please store batteries and diagnostic scale at a place children cannot reach. In case a battery has been swallowed, please consult a physician immediately.

Keep children away from the packaging materials - danger of asphyxiation!

Do not recharge the batteries. Do not take them to pieces. Do not throw them into a fire. Do not short-circuit them.

III. Setup & Operating functions

Changing the batteries

- Please open the cover of the battery case at the back of the device.
- Insert 4 AAA batteries (1,5 V). Please pay attention to the polarity.
- Close the cover of the battery case.

Note: In case the scale does not operate please remove the batteries and insert them again.

Your **smartLAB® fit** body analyser scale has a symbol indicating that you have to change the batteries. In case the batteries are almost empty, the symbol “Lo” can be seen on the display. Please change the batteries.

Please remove the batteries in case you are not going to use the scale for one month or more in order to avoid an automatic discharge of the batteries.

The scale as well as the batteries must be disposed of according to the local regulations.

Beginning of operation / „quick weighing “

- Use the **smartLAB® fit** body analyser scale bare-footed. Do not move so that you can stand on both legs. The scale starts measuring immediately and finally shows your body weight on the display. When you leave the tread of the scale it deactivates automatically.

Setting up of the user data

- In order to measure the adipose content and other parameter you first have to insert your individual user data. You can store up to 8 user profiles and recall the measured data.

Please press “SET” (on the side of your scale) to create a new user profile.

Note: The weight scale can be set by the touch function, which can be found on the right top side of the scale. The touch function represents the buttons, which are abreast of the Set writing and the arrows at the side of the scale.

User profile:

The display now shows several illuminated preadjustments. The symbol “P1” appears on the display (bottom, left) for user profile 1 (see picture 1). You can confirm this choice by pressing “SET” or you can create a new user profile by increasing the number (arrow keys at the side of the scale).

Pic. 1

Gender:

Press „SET“ to confirm your choice and then choose your gender. The display shows the symbol for a man. If you want to change the gender please use one of the arrow keys until you can see the symbol for a woman. Press “SET” to confirm your choice.

Body height:


Please insert your body height by using the arrow keys. You can insert body heights from 100 cm up to 250 cm. When you continue to press the buttons the indicated values go faster. Please confirm your value by pressing “SET”.

Age:

The display shows the value for your age. Please use again the arrow keys to insert your age and confirm by pressing “SET”.

Activity index:

Then the symbol for the choice of the activity index appears in the left display segment (see picture 2).

Pic. 2

The value defines a certain degree of activity that is important for a medium-term or a long-term consideration. You can choose between 4 activity degrees for each profile by using the arrow keys. Please confirm your choice and press "SET".

Defining the different degrees of activity:

- Degree 1: No physical exercise
- Degree 2: Only little exercise such as promenades, working in the garden or gymnastics.
- Degree 3: Medium physical exercise (30 minutes physical exercise from 2 up to 4 times per week minimum)
- Degree 4: High physical exercise (30 minutes physical exercise from 4 up to 6 times per week minimum)

After having finished the set up of your data the symbol "0,0" appears on the display. The **smartLAB[®] fit** body analyser scale is ready to start the measurement.

Measuring

After having determined all parameter and user profiles body weight, adipose and the other values can be determined.

- Please press „SET“ and choose your user profile by using the arrow keys. Do NOT confirm by pressing "SET" since you would have to insert and confirm again all parameters such as gender, body weight etc.
- The measurement starts automatically after a few seconds. First of all your **smartLAB[®] fit** body analyser scale detects your body weight (left display segment). After that the measurement of the other parameter will start. You will see three circles on the display segments concerned (please refer to chapter display & functions in this user manual). The scale is detecting the other parameters – content of water, content of muscles and bone mass (see picture 3).

Pic. 3

- Finally all values are displayed (see picture 4). After a few seconds the indication of the body weight changes into the indication of the calorie consumption (see picture 5), then into the indication of the activity (in kcal). In addition the activity index is

indicated by the symbol “activity index” (see picture 6). At the end your **smartLAB® fit** body analyser scale displays again the parameters for your user profile which are the basic data for your measured values.

Pic. 4



Pic. 5

Pic. 6

Optional: Weighing on a carpet

You can use the **smartLAB® fit** with the extension feet on a carpet. For this please follow these steps:

1. Remove the non-slip pads from the scales feet
2. Fit the extension feet firmly onto the scales feet. Now you can use the scale on a carpet.

Interpreting the results

The following guiding values help you to interpret the measured values of your **smartLAB® fit** body analyser scale. These values are only an approach you should eventually discuss and interpret with your physician.

Adipose content: ●●●

The following values are guiding values concerning your adipose content. Please contact your physician for further information.

Age	Woman				Man			
	undernourished	healthy / normal	overweigh	too much adipose	undernourished	healthy / normal	overweigh	too much adipose
10 - 12	< 12	12 - 22	22 - 30	> 30	< 8	8 - 18	18 - 24	> 24
13 - 18	< 15	15 - 25	25 - 33	> 33	< 8	8 - 18	18 - 24	> 24
19 - 30	< 20	20 - 29	29 - 36	> 36	< 8	8 - 18	18 - 24	> 24
31 - 40	< 22	22 - 31	31 - 38	> 38	< 11	11 - 20	20 - 26	> 26
41 - 50	< 24	24 - 33	33 - 40	> 40	< 13	13 - 22	22 - 28	> 28
51 - 60	< 26	26 - 35	35 - 42	> 42	< 15	15 - 24	24 - 30	> 30
60+	< 28	28 - 37	37 - 47	> 47	< 17	17 - 25	26 - 34	> 34

Note: Athletes often have lower values than those of the table above. The values can be underneath the guiding values depending on the sport, the intensity of exercise and the physical condition. But please note that extremely low values can be an indication for health hazard.

Content of water of your body:

The content of water inside your body is usually within the following ranges (only the normal range is indicated):

Age	Woman	Man
10 - 12	> 60	> 64
13 - 18	> 58.5	> 63.5
19 - 30	> 56	> 62.5
31 - 40	> 53	> 61
41 - 50	> 52	> 60
51 - 60	> 51	> 59
60+	> 50	> 58

Note: Adipose has relatively few water. For this reason the content of water in case of a high content of adipose can be below the guiding values above. Regarding endurance athletes the values can be higher than the guiding values because of the low adipose and the high content of muscles.

Content of muscles:

The content of the muscles usually stays within the following ranges (only the normal range is indicated):

Age	Woman	Man
10 - 99	> 34	> 40

Bone mass:

Our bones as well as our whole body are subject to the ageing process and all its consequences. During the childhood the bone mass is increasing fast and reaches its maximum at the age of 30-40. Then the bone mass is decreasing. You can influence the ageing process by a healthy nutrition (esp. calcium and vitamin D) and by a continuous physical exercise. In addition to that you can sustain your bones by a specific muscle formation.

 Bone mass and bone density are different factors. Please do not confuse them. The latter can only be detected by a medical examination (e. g. computer tomography, ultra sound). For this reason the scale cannot be used to identify changes or the bones or their durability.

The bone mass can hardly be influenced but is changing according to factors like weight, body height, age or gender.

	Woman			Man		
Body weight	< 45 kg	45 - 60 kg	> 60 kg	< 60 kg	60 - 75 kg	> 75 kg
Bone mass	1,8 kg	2,2 kg	2,5 kg	2,5 kg	2,9 kg	3,2 kg

Basal Metabolic Rate (BMR):

BMR means the amount of energy the body needs for its basal vital functions when the body is in the state of absolute calmness (e. g. when you stay in bed for 24 hours). This value mainly depends on factors like weight, body size and age. This value is displayed on your **smartLAB® fit** body analyser scale. The unit is kcal/Tag (kcal/day). The basis for the calculation is the science-based Harris-Benedict-Formula.

This is the amount of energy your body needs at least in the form of nutrition. If your body gets less energy for a longer time this can affect your health negatively.

Active Metabolic Rate (AMR):

The active metabolic rate means the amount of energy the body needs daily in the state of activity. The energy consumption increases with an increasing physical exercise. Your **smartLAB® fit** body analyser scale has 4 different activity indices (refer to chapter “Setting up of the user data”) to determine the appropriate degree of activity.

In order to maintain the current weight the consumed energy must be replaced in the form of food and drink. In case this energy will not be replaced the body uses the stored fat – the body loses weight. In case the body gets more energy than the required AMR the body is not able to use this energy – it will be stored inside the body and you gain weight.

Time as a factor:

Please note that only a long-term result counts. Short-term changes of your weight within a few days are absolutely normal and are often due to changes of the content of water. Only middle-term or long-term changes are related to the content of your adipose or your muscles. In case your weight decreases for a short time and the content of the adipose of your body increases or stays the same you only lost water (e.g. after a physical exercise or after having visited a sauna). If your weight increases in the medium term and your content of adipose decreases or stays the same you might have build up more muscle mass.

Note: Adipose, content of water or muscle mass may not be added because the muscle tissue also contains components of water.

IV. Miscellaneous

Maintenance

- Please clean your device regularly by a damp cloth. Do not use caustic cleaning solutions. Do not put the device under fluent water.
- Protect the device from impacts, humidity, dust, chemicals or fluctuations of temperature. Do not store near heat sources such as ovens or heating installations.
- NEVER take the scale to pieces by yourself. In this case the warranty will not be valid anymore. Only the manufacturer or the authorized trader is allowed to repair the scale. Please check the batteries before any complaint and change them if necessary.
- The **smartLAB® fit** diagnostic scale is in accordance with the EC directive 2004/108. Please feel free to contact the manufacturer or the customer service for further questions regarding the applications of the scale.

Regulation and Safety Notice

WEEE note

The WEEE (Waste Electrical and Electronic Equipment) Directive, that came into effect as European law on 13th February 2003, led to a major change in the disposal of electrical equipment.

The primarily purpose of this directive is to prevent electronic waste (WEEE), while encouraging for reuse, recycling and other forms of reprocessing to reduce waste.

The WEEE logo on the product and the packaging indicates that it is not allowed to dispose of the product in the ordinary household waste. It is your responsibility to deliver all the disused electrical and electronic equipment to the respective collection points. A separate collection and proper recycling of electronic waste helps dealing economically with natural resources. Furthermore, the recycling of electrical waste is a contribution to the conservation of our environment and thus the health of all people. More information concerning the disposal of electrical and electronic devices, reprocessing and the collection points you can get from the local authorities, waste disposal companies, from retailers and manufacturers of the device.

RoHS compliance

This product complies with Directive 2002/95/EC of the European Parliament and the Council of 27th January 2003 regarding the limited use of dangerous substances in electrical and electronic equipment (RoHS) and its variations.

Warranty

HMM Diagnostics GmbH puts emphasis on a high quality standard. For this reason HMM Diagnostics GmbH gives a warranty of 2 years. You can extend your warranty from 3 up to 5 years by enrolling your product at HMM Diagnostics GmbH. Please register online under following website: www.hmm.info/en/registration.

Parts subject to regular wear and tear, batteries etc. are not subject to the warranty.

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